

# Dual Water Washout Tester 19650-0

## Water Washout

ASTM D1264; IP 215; DIN 51807-2; ISO 11009

- Twin bearing and housings for duplicate or independent testing
- Two separate water tanks and jets
- All parameters are independently adjustable
- Testing at 38 and 79 °C (100 and 75 °F) as per ASTM D1264
- Magnetic bearing drive for easy disassembly and removal
- Automatic test procedure and display of test parameters
- Spacious and easy to clean water tank
- All water wet parts are resilient to salt water



Technical Specifications	
Operation	
Rotational speed	600 and 3780 rpm (variable options available)
Temperature range	Ambient to 80 °C
Power Requirements	
Voltage	220/240V 50Hz
Physical	
Dimensions (HxWxD)	600 x 850 x 520 mm
Weight	55 kg

• Greases •